

BriteView Color Code Laryngoscope Blades

Choose the BriteView™ color-code disposable blade system and know you have the right tool for the job

- **Super bright light** for maximum visibility compared to bent clear plastic rod design

- **Stainless Steel from tip to heel** prevents flexing or breaking under stress



- **Single patient use** prevents cross contamination and eliminates cleaning costs
- **Color Coded** for instant size recognition
- **Fits all standard handles** which are a fraction of the cost of fiber-optic types

BriteView™ Color Code Blade System

Designed by EMS professionals with real-life street experience, the BriteView™ Laryngoscope Blade Series provides a strong, bright, laryngoscope blade, color coded for instant size recognition. Disposable stainless steel construction means you don't have to compromise strength and reliability for the cleaning and maintenance-free benefits of a single patient use blade.



The unique color coded heels make it a snap to select the correct size blade for the job at hand. The blade sizes are also etched on the side of the blade in a large format for easy selection when laid out on an instrument tray. The innovative package graphics are large and clear for easy identification in an emergency airway bag, under low light conditions or when restocking your vehicle's inventory.

If you want a disposable blade that won't flex, won't break off the handle and will provide a bright, reliable light source, then select the BriteView blade series today. Specify the BriteView Color-Code blade on your next bid request. Contact your authorized Hartwell Medical dealer or give us a call at **1-800-633-5900**.

Available in Macintosh and Miller sizes



Model

LB 42MIL0-DS
 LB 42MIL1-DS
 LB 42MIL2-DS
 LB 42MIL3-DS
 LB 42MIL4-DS
 LB 42MAC1-DS
 LB 42MAC2-DS
 LB 42MAC3-DS
 LB 42MAC4-DS

Description

Miller 0 - White
 Miller 1 - Pink
 Miller 2 - Yellow
 Miller 3 - Orange
 Miller 4 - Lime Green
 Macintosh 1 - Pink
 Macintosh 2 - Yellow
 Macintosh 3 - Orange
 Macintosh 4 - Lime Green

HARTWELL MEDICAL

6352 Corte del Abeto, Suite J
 Carlsbad, CA 92011-1408
 800-633-5900
 760-438-5500 Tel
 760-438-2783 Fax
www.HartwellMedical.com

Creating Quality Products for Life