

CATCH - ALL™ - The Versatile Equipment Storage Tray



- Simple and Easy to Use
- Strong, Durable Fabric
- Waterproof Material Will Not Absorb Blood or Body Fluids
- Unique Packaging Pouch Becomes a See-Thru Storage Pouch
- Specific Models Style Cots for Roll-In and Lift-In Style Cots
- Simple Soap and Water Clean-Up

The CATCH-ALL™ solves the problem of where to keep essential equipment like gloves, stethoscope, BioHoop® Bags and a pocket mask. It's also great for patient belongings, medicine vials, X-Rays and your run report clipboard.

The Hartwell Medical CATCH-ALL™ provides an easy and convenient storage area for small medical equipment that is frequently used on a majority of EMS calls. The flexible, vinyl-coated material is lightweight, yet durable so it doesn't add unnecessary weight to your cot. The CATCH-ALL™ attaches to the cot's main frame underneath the backrest in less than two minutes. Four models of the CATCH-ALL™ have been created to fit the majority of ambulance cots that

are currently used in the field today. Whether you need storage for a pocket mask and stethoscope, or a handy shelf to transport the patient's clothing or other personal belongings, the CATCH-ALL™ meets the needs of today's EMS providers. Remember to specify the CATCH-ALL™ model number when ordering.

The CATCH-ALL™ is designed to fit your cot, whether it is a lift-in style cot or a roll-in style cot. Make sure to specify which cot you are currently using to ensure that your CATCH-ALL™ will fit each and every cot in your fleet. Models CA 3000 and CA 4000 are shown above.



Ferno is a registered trademark of Ferno-Washington, Inc. Stryker and MX-PRO are registered trademarks of Stryker Corporation

ORDERING INFORMATION:

COT MODEL	CATCH-ALL MODEL	COLORS
Ferno® Model 30	CA 2000	Black
Ferno® Model 35A	CA 3000	Black
Ferno® Model 29M	CA 4000	Black
Ferno® Model 93ES	CA 4000	Black
Stryker® MX-PRO®	CA 6000	Black

SPECIFICATIONS

DIMENSIONS:

CA 2000:	20.5" x 21.5"
CA 3000:	18.0" x 20.0"
CA 4000:	18.0" x 20.0"
CA 6000:	20.0" x 21.5"

Weight: 5 to 7 ounces (Depending on model)

Materials:

Fabric:	Vinyl/Polyester
Fasteners:	Nylon

HARTWELL MEDICAL

6354 Corte del Abeto, Suite F

Carlsbad, CA 92011-1479

1-800-633-5900

(760) 438-5500 Phone

(760) 438-2783 Fax

www.HartwellMedical.com

Creating Quality Products for Life

NEW!
CA6000 for
Stryker® MX-PRO®