

# BriteGrip™

**LIGHTWEIGHT DISPOSABLE**

## Laryngoscope Handle

**The performance and power of a reusable laryngoscope handle without the weight, contamination hazard or hassle**

Weighing in at just 3 ounces, the Hartwell Medical, single-use disposable laryngoscope handle provides the performance and power of a reusable at a fraction of the weight. This innovative handle offers numerous advantages to substantially reduce the potential for cross-contamination known to exist with reusable devices. Plus, you will save valuable time and hassle never needing to sterilize and inspect your handles again. Eliminates the cost and time spent on battery inspections, maintenance and inventory.

The sturdy construction of the BriteGrip insures reliability in the most adverse field conditions. An ideal companion product to BRITEVIEW® and GRANDVIEW® laryngoscope blades.

### Features and Benefits

- **Incredibly Lightweight, Only 3 oz.**
- **Batteries Included, Ready-To-Go!**
- **Single-Use Design & Disposable**
- **Reduces Potential Cross-Contamination**
- **Sturdy Construction**
- **Powerful, Reliable And Recyclable**
- **Eliminates Equipment Sterilization**
- **Compatible With Standard Blades That Conform To ASTM (F965) And ISO (7376-1)**

Weighs Only 3 oz.

*BriteGrip Part Number LB 1002-DS*

