

# CATCH-ALL™ The Versatile Equipment Storage Tray




- Simple and Easy to Use
- Strong, Durable Fabric
- Waterproof Material Will Not Absorb Blood or Body Fluids
- Unique Packaging Pouch Becomes a See-Thru Storage Pouch
- Specific Models for Roll-In and Lift-In Style Cots
- Simple Soap and Water Clean-Up

**The CATCH-ALL™ solves the problem of where to keep essential equipment like gloves, stethoscope, BioHoop® Bags and a pocket mask. It's also great for patient belongings, medicine vials, X-Rays and your run report clipboard.**

The Hartwell Medical CATCH-ALL provides an easy and convenient storage area for small medical equipment that is frequently used on a majority of EMS calls. The flexible vinyl-coated material is lightweight, yet durable so it doesn't add unnecessary weight to your cot. The CATCH-ALL attaches to the cot's main frame underneath the backrest in less than two minutes. Various models of the of the CATCH-ALL have been created to fit many of the ambulance cots that are currently used in the field today. Whether you need storage for a pocket mask and stethoscope, or a handy shelf to transport the patient's clothing or other personal belongings, the CATCH-ALL meets the needs of today's EMS providers. Remember to specify the CATCH-ALL model number when ordering.

The CATCH-ALL is designed to fit your cot, whether it is a lift-in style cot or a roll-in style cot. Make sure to specify which cot you are currently using to ensure that your CATCH-ALL will fit each and every cot in your fleet. Model CA 3000 and CA 4000 are shown here.



## ORDERING INFORMATION

COT MODEL	CATCH-ALL MODEL	COLOR
Ferno® Model 35A	CA 3000	Black
Ferno® Proflex Series	CA 3600	Black
Ferno® Model 29M	CA 4000	Black
Ferno® Model 93ES	CA 4000	Black
Stryker® MX-Pro®	CA 6000	Black

## SPECIFICATIONS

### DIMENSIONS

CA 3000	18.0" x 20.0"
CA 3600	18.0" x 20.0"
CA 4000	18.0" x 20.0"
CA 6000	20.0" x 21.5"

Weight: 5 to 7 ounces (Depending on model)

Materials:

Fabric: Vinyl/Polyester

Fasteners: Nylon

Ferno is a registered trademark of Ferno-Washington, Inc. Stryker and MX-PRO are registered trademarks of Stryker Corporation