

# BriteView®

**COLOR-CODED DISPOSABLE**

## Laryngoscope Blade System

**Designed by EMS professionals with real-life street experience. The BriteView® Disposable Laryngoscope blade provides a no flex blade, a bright light and is color-coded for instant size recognition.**

Stainless steel construction means you do not have to compromise strength and reliability for the cleaning and maintenance-free benefits of a single patient use laryngoscope blade.

The unique color-coded heel makes it a snap to select the correct size blade for the job at hand. The blade sizes are also etched on the side of the blade in a large format for easy selection when placed on an instrument tray. The innovative package graphics are large and clear for easy identification in an emergency airway bag, under low-light condition or when restocking your vehicle's inventory.

If you want a disposable blade that will not flex, will not break off the handle and will provide a bright reliable light source, then select the BriteView Blade System.



*Available in Miller and Macintosh sizes*

### Features and Benefits

- **Stainless Steel From Tip To Heel Prevents Flexing Or Breaking Under Stress**
- **Color-Coded For Instant Size Recognition**
- **Single Patient Use Prevents Cross Contamination**
- **Eliminates Cleaning Costs**
- **Super Bright Light Gives Maximum Visibility**
- **Tear-Pack Is Easy To Open With Gloves On**
- **Fits All Standard Handles**
- **Available In Miller And Macintosh Styles**

| MODEL        | DESCRIPTION              |
|--------------|--------------------------|
| LB 42MIL0-DS | Miller 0 - White         |
| LB 42MIL1-DS | Miller 1 - Pink          |
| LB 42MIL2-DS | Miller 2 - Yellow        |
| LB 42MIL3-DS | Miller 3 - Orange        |
| LB 42MIL4-DS | Miller 4 - Lime Green    |
| MODEL        | DESCRIPTION              |
| LB 42MAC1-DS | Macintosh 1 - Pink       |
| LB 42MAC2-DS | Macintosh 2 - Yellow     |
| LB 42MAC3-DS | Macintosh 3 - Orange     |
| LB 42MAC4-DS | Macintosh 4 - Lime Green |