



# THE ULTIMATE Collection Bag

**The BioHoop® is a versatile and leak-proof bag that solves your waste containment needs and much more**

Featuring the largest opening of any compact collection bag, the BioHoop makes your job easier and reduces your clean-up time. Semi-rigid rim keeps the bag open, providing a large target for the patient that will not collapse. Large 8" x 13-3/4" bag holds up to 3 liters. Long, protective cuff on both sides protects your hands from fluids. Cinchable strap allows instant closure with a quick pull, eliminating spills and odor. Unique hook allows the bag to hang within easy reach on ambulance cot rails, hospital beds or wheelchairs.

Clear front with graduated markings makes it easy for hospital personnel to examine contents. Tamper resistant bag is perfect for emesis, urine, blood spills, soiled dressing, patient belongings like jewelry and cell phones, evidence or other waste material.



*BioHoop with Hook Part Number BH 1100H  
BioHoop without Hook Part Number BH 1100*



Emesis



Dressings & Scene  
Garbage



Patient Belongings  
& Evidence



Medications



Urine



*Easily and safely manage scene garbage or patient belongings quickly and efficiently.*



## INNOVATIVE DESIGN

Semi-rigid rim keeps the bag open, then closes instantly with a zip of the strap.



## QUALITY MATERIALS

Long, protective cuff at the rim protects your hands from body fluids.



## MORE THAN JUST A PLASTIC BAG

The patient's personal items can be secured in the tamper-resistant bag and attached to their clothes with the unique hook clip to ensure nothing is lost.



## SMART FEATURES

Clear front makes it easy to examine the contents and helps comply with OSHA 29 CFR Part 1910.1030.