

CombiCarrierII[®]

Head Immobilizer

Unique Design

Designed specifically for use with the award-winning features of the CombiCarrierII[®], this head immobilizer allows for easy operation of the latch mechanism on the CombiCarrierII during application and removal from the patient.

Saves You Time

You can leave the low-profile base plate on the CombiCarrierII when storing it and save valuable on-scene time during an emergency call.

Security and Easy Patient Access

Compact, anatomically shaped head supports minimize lateral movement and maximize access to the patient's neck and ears without compromising stabilization.

Durable and Versatile

Constructed of vinyl-dipped NBR, the head immobilizer is X-ray compatible, easily cleaned, and requires minimal storage space. It also fits most other 16" backboards in use today and is ideal for water rescues.

Reliable Quality

Hartwell Medical is known for its quality workmanship, outstanding customer service and warranty program. Should you have any questions about the use or care of any of our products, our customer service personnel will be happy to help you.

Use the model numbers below to order your head immobilizer or replacement parts.

<u>MODEL #</u>	<u>DESCRIPTION</u>
CC 3400H	CombiCarrierII Head Immobilizer
<u>REPLACEMENT PARTS</u>	<u>DESCRIPTION</u>
CC 3400-9-100	CombiCarrierII Head Immobilizer Base Plate (each)
CC 3400-9-200	CombiCarrierII Head Immobilizer Head Support Block (pair)
CC 3400-9-300	CombiCarrierII Head Immobilizer Head Strap (each)



HARTWELL MEDICAL

800-633-5900 • 760-438-5500 • 760-438-2783 (Fax)

www.HartwellMedical.com

Creating Quality Products for Life

